

Vodivostní standardy ROTI® CALIPURE

Standardní roztoky pro kalibraci konduktometrů

- Dodávány s certifikátem šarže
- Traceable to SRM from NIST

Technické parametry

Přesnost: $\pm 1 \%$ (25 °C)



Standard při 25 °C (μS)	Balení	Katalogové číslo
5	500 ml	R.2462.1
10	250 ml	R.0810.1
10	500 ml	R.0810.2
20	500 ml	R.0811.1
50	500 ml	R.0812.1
84	500 ml	R.2432.1
84	1 l	R.2432.2
100	500 ml	R.2463.1
147	500 ml	R.2433.1
147	1 l	R.2433.2
500	500 ml	R.2464.1
1 000	500 ml	R.2465.1
1 413	500 ml	R.2434.1
1 413	1 l	R.2434.2
5 000	500 ml	R.0814.1
10 000	500 ml	R.2466.1
12 880	500 ml	R.2435.1
12 880	1 l	R.2435.2
20 000	500 ml	R.0822.1

50 000	500 ml	R.0831.1
100 000	500 ml	R.0832.1
500 000	500 ml	R.2467.1